



SOLIDEX[®]-ISOEx Human CD3 Isolation Nanobeads (Column-Based)

Product Instruction Manual

Cat. No.	GM-Tg-hg-T87075-iso-nanoIMB
Product Name (Column-Based)	SOLIDEX [®] -ISOEx Human CD3 Isolation Nanobeads
Storage Temperature	2-8 °C (Do not freeze)

Technical Parameters

Classification	Description
Cell type	CD3+ Pan-T Cell
Isolation method	Positive selection
Labeling type	Direct
Magnetic nanobeads type	Non-releasable
Format	Column-based

Product Introduction

T cells are a critically important class of immune cells in the human body, accounting for 45% to 70% of peripheral blood mononuclear cells (PBMCs). They play a vital role in anti-infection responses, anti-tumor immunity, and immunomodulation. Depending on the specific application scenario, we support three isolation strategies: positive selection, negative selection, and depletion.

By utilizing SOLIDEX[®]-ISOEx Human CD3 Isolation Nanobeads (Column-Based), column-based isolation achieves superior cell purity. Through the specific capture facilitated by the magnetic matrix within the column, high-purity target cells can be effectively obtained. This method is highly compatible with automated workflows and facilitates various downstream applications, such as investigating T cell cytotoxicity, activation processes, HIV infectivity, signal transduction mechanisms, and surface marker expression.

Product Components and Specifications

Component	Catalog No.	10 Test (1×10 ⁸ cells)	50 Tests (5×10 ⁸ cells)	100 Tests (1×10 ⁹ cells)
SOLIDEX [®] -ISOEx Human CD3 Isolation Nanobeads (Column-Based)	GM-Tg-hg-T87075-iso-nanoIMB	100 µL	0.5 mL	1.0 mL

Note: This product is for research and development use only.

Shelf Life: Store at 2-8 °C, protected from light, and do not freeze. Under these conditions, the product is valid for 6 months.

Reagents and Equipment Required

A. Cell isolation Column

For positive selection or depletion, recommended for use with GeneMedi SOLIDEX[®]-ISOEx cell isolation columns: M Column for standard throughput, Cat. No.: **GMP-ISOEx-Column-M**; L Column for high throughput, Cat. No.: **GMP-ISOEx-Column-L**. Comparable columns from other mainstream brands are also compatible.

B. Magnetic Separator

C. Cell isolation buffer

Phosphate-buffered saline (PBS), pH 7.2, containing 0.5% bovine serum albumin (BSA) and 2 mM EDTA. (user-supplied).

Note:

- (1) BSA can be replaced by other proteins such as human serum albumin (HSA), human serum, or fetal bovine serum (FBS).
- (2) Degas the Cell isolation buffer before use, as air bubbles may block the Column.
- (3) Keep the Cell isolation buffer cold (2-8°C).

Protocol

A. Cell Labeling

- a. Count the peripheral blood mononuclear cells (PBMCs). Add 50 μ L of cell isolation buffer per 1×10^7 cells and resuspend the cells.
- b. Add 10 μ L of SOLIDEX®-ISOEx Human CD3 Isolation Nanobeads per 1×10^7 cells. Mix gently but thoroughly, and incubate at 4°C for 15 minutes.
- c. Add 500 μ L of cell isolation buffer per 1×10^7 cells. Centrifuge and discard the supernatant. Resuspend the cell pellet in 3 mL of cell isolation buffer. The cells are now ready for the subsequent cell isolation step.

Note:

- a. The nanobeads must be thoroughly mixed prior to use by pipetting up and down.
- b. The cell isolation buffer must be pre-cooled to 2-8°C or on ice.

B. Cell Isolation

- a. Place the LS isolation column on the magnetic separator and rinse with 3 mL of cell isolation buffer.
- b. Once the cell isolation buffer has completely drained, apply the 3 mL cell suspension to the column. Place a clean 15 mL centrifuge tube underneath to collect the flow-through.
- c. After the liquid has drained from the column, wash the column once by adding 5 mL of cell isolation buffer. Then, remove the isolation column from the magnetic separator and place it into a clean 15 mL centrifuge tube. Add 5 mL of cell isolation buffer and use the plunger to directly flush the liquid out of the column; this collected fraction contains the CD3+ T cells.

Notes

- A. Avoid freezing during use and storage of the Nanobeads.
- B. It is recommended to use low-binding pipette tips and centrifuge tubes to prevent loss of Nanobeads due to adsorption.
- C. Before aspirating the Nanobeads, mix them gently. Avoid bubble formation during mixing.
- D. This product is used for research use only.

Validation Data from GeneMedi




Figure 1. Isolation of high-purity T cells achieved by the SOLIDEX®-ISOEx Human CD3 Isolation Nanobeads (Column-based). To evaluate the purity of the isolated T cells, CD3⁺ T cells were isolated from human peripheral blood mononuclear cells (PBMCs). Cells pre- and post-isolation were labeled with a CD3-APC antibody and PI (Propidium Iodide) for flow cytometric analysis. The purity of CD3⁺ T cells pre-isolation and post-isolation was 52.6% and 95.0%, respectively.