



PurProX® AAVEasy

Rapid-to-purify Spin Column

『Easy』 Your AAV Purification



Easy Spin-Purification
Ultracentrifugation-FREE



Serotype All-in-One
One solution for Most AAV serotypes



High Recovery Rate
>90% AAV recovery



Easy-Handling
within 30 min



High throughput
Dozens of samples purified at once



1000xGMXbuffer01 Lower plug Resin tube Resin tube + Lower plug Sample tube

PurProX® AAVEasy + AAVFull = High-Purity, High Full-Capsid AAV



Scan →



Brochure Download



Sample request

PurProX® AAVEasy rapid-to-purify Spin Column — Efficient Serotype-Specific

Cat No,	Product Name	Serotype
GMV-PurProX-AAVEasy-U-1	PurProX® AAVEasy rapid-to-purify Spin Column-1	Versatile resin for high-capacity, multi-serotype binding (e.g., AAV2/5/6/8/DJ/7m8)
GMV-PurProX-AAVEasy-U-2	PurProX® AAVEasy rapid-to-purify Spin Column-2	Optimized for AAV9 and variants like AAV-PHP series

PurProX® AAVFull Enrichment Kit — Specifically clears Empty AAV Capsids




Cat No,	Product Name	Feature
GMV-PurProX-AAVFull-U-1	PurProX® AAVFull Enrichment Kit	High Binding Capacity: $2-5 \times 10^{13}$ VP/mL



PurProX[®] AAV Purification Solution

PurProX[®] AAVEasy purification process



PurProX[®] AAVEasy Spin Column-1 for AAV Serotype 2, 5, 6, 8, DJ, and 7m8 Purification

PurProX[®] AAVEasy Spin Column-2 Designed for AAV9 Purification

